



Digital Thermometer

GAMMEX TM-99A

Minor shifts in developer temperature can have a detrimental effect on the film contrast and density. In order to achieve and maintain appropriate film speed, film contrast and film fog levels, the developer temperature must be monitored on a regular basis. To maintain this temperature a high quality, accurate thermometer must be used, such as the Gammex TM-99A.

This thermometer, with its fast acting probe, degrees Celsius or Fahrenheit and low battery indication is the perfect choice for all processor quality control procedures.

SPECIFICATIONS

- Range -40 to 150° C (-40 to 302° F)
- Accuracy ±0.2° F from -40 to 100° F;
±2% from 100 to 300° F
- Power 9 V Alkaline Battery (included)
- Size 3x7.5x18 cm (1.2x3x7 in)
- Weight 255 g (9 oz)

