



Routine Mammographic QCK it

GAMMEX 183

The Gammex1 83 Routine Mammographic QCK it contains the tools, instructions, and recording forms needed for a film-screen mammography program that meets MQSA requirements. The instruments evaluate image quality, compression force, film-screen contact, processor performance, and film hyperfotometry. In addition, the kit includes a camel hair brush and lint-free cloths for maintaining mammographic screens.

SPECIFICATIONS

Gammex1 83 kit contains:

- Densitometer MA0025
- Sensitometer M A5034
- Digital Thermometer T M-99A
- Mammographic Accreditation Phantom 156
- 4mm acrylic Disk 158
- Mammographic Film/Screen Test Tool 157A
- Breast Compression Device 163
- Fixer Retention kit 166B
- Lint Free Cloths 167
- Camel Hair Brush 168
- ACRM mammographic Quality Control Manual
- Foam-lined Case 081A

Size..... 24x42x63.5cm (9.5x16.5x25.1 in)
 Weight..... 9kg(20 lbs)

WWW.GAMMEX.COM