



SPECIFICATIONS

MODEL 136500 INCLUDES

Qty	Description
1	Head Vertex
1	Head & Shoulders Anterior part
1	Head & Shoulders Posterior part
1	End plate, ABS black
4	Solid Rod w/BB @ ISO center, 330 mm long
3	Sleeve CV511C, 115 mm long
1	Plug PL-CV511C solid, 330 mm long
1	Plug PL-CV511C solid w/BB, 330 mm long
1	Plug PL-CV511C solid to fit Cavity Rod @ ISO center, 217 mm long
1	Plug PL-CV511C solid to fit Cavity Rod @ ISO center + 50 mm, 167 mm long
1	Spacer Plugs CV511C Kit, total length 370 mm x 7mm Ø
1	Spacer Plugs CV511C Kit, total length 130 mm x 13mm Ø
1	Cavity Rod CV511C for Ion Chamber @ ISO location, 320 mm

Qty	Description
1	Cavity Rod CV511C for Ion Chamber @ ISO+50 mm location, 320 mm
5	Water vial
1	Screw driver
12	Nylon screws ¼-20x1-1/4"
6	Black nylon thumbscrews ¼-20x1/2"
8	Nylon pins 1/4x1"
1	Black nylon strap with a buckle
1	User's Guide (electronic copy)
1	Carry case, foam lined

Dimensions: 36 cm (W) x 36 cm (L) x 22cm (H)

Weight: 10.2 kg without rods

Materials: CIRS proprietary epoxy resins.

©2016 Computerized Imaging Reference Systems, Inc. All rights reserved. All brand names, product names or trademarks belong to their respective holders. Specifications subject to change without notice. Publication: 1365 DS 120418



Computerized Imaging Reference Systems, Inc. has been certified by UL DQS Inc. to (ISO) 13485:2003. Certificate Registration No.10000905-MP23..

PEO B.V.

info@gotopeo.com
www.gotopeo.com

The Netherlands

Havenweg 16, 6603 AS Wijchen
+31 (0)24 648 86 88

Belgium

Watermolenstraat 2, B-2910 Essen
+32 (0)3 309 32 09

CoC 34107894

VAT NL807859151B01

IBAN NL29 RABO 0356 1960 46

BIC RABONL2U

